

We have added Passivation to our process capabilities.

AMS 2700 or QQ-P-35

Method 1 = Nitric

Method 2 = Citric

Capabilities:

Method 1, types 1,4,6

Method 2 (All processed the same)

Testing by Method 102 Copper Sulfate only

ASTM A967

Capabilities:

Types Nitric 2, Citric 3, Nitric 5

Testing by Practice D, F, G

If you would like to discuss this process, please contact me or your
Customer Service Representative.

Regards,

Cliff Maddox

Cliff Maddox
Quality Manager
602-276-1231 x5